

Product Datasheet

Rat Ab1-40 ELISA Kit (Ready to Use) (orb554617)

Description	Rat Amyloid Beta Peptide 1-40 (Ab1-40) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	11.6pg/mL
Expiration Date	Please enquire.