

Product Datasheet

Rat APBB3 ELISA Kit (Ready to Use) (orb554620)

Description	Rat Amyloid Beta Precursor Protein Binding Protein B3 (APBB3) Sandwich ELISA KitReady to Use: Pre-di...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8R1C9
Sensitivity	0.062ng/mL
Entrez	225372
Expiration Date	Please enquire.