

Product Datasheet

Rat AGGF1 ELISA Kit (Ready to Use) (orb554623)

Description	Rat Angiogenic Factor With G Patch And FHA Domains 1 (AGGF1) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Rat
Range	0.468-30ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F1MAK0
Sensitivity	0.19ng/mL
Entrez	310005
Expiration Date	Please enquire.