

Product Datasheet

Rat AGTR2 ELISA Kit (Ready to Use) (orb554638)

Description	Rat Angiotensin II Receptor 2 (AGTR2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P35351
Sensitivity	0.051ng/mL
Entrez	24182
Expiration Date	Please enquire.