

Product Datasheet

Rat APOA2 ELISA Kit (Ready to Use) (orb554650)

Description	Rat Apolipoprotein A2 (APOA2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for con...
Reactivity	Rat
Range	6.25-400ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04638
Sensitivity	2.57ng/mL
Entrez	25649
Expiration Date	Please enquire.