

Product Datasheet

Rat ARNTL ELISA Kit (Ready to Use) (orb554675)

Description	Rat Aryl Hydrocarbon Receptor Nuclear Translocator Like Protein (ARNTL) Sandwich ELISA KitReady to U...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9EPW1
Sensitivity	0.057ng/mL
Entrez	29657
Expiration Date	Please enquire.