

Product Datasheet

Rat ATP1a1 ELISA Kit (Ready to Use) (orb554684)

Description	Rat ATPase, Na ⁺ /K ⁺ Transporting Alpha 1 Polypeptide (ATP1a1) Sandwich ELISA KitReady to Use: Pre-dil...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P06685
Sensitivity	0.062ng/mL
Entrez	24211
Expiration Date	Please enquire.