

Product Datasheet

Rat ATP1b4 ELISA Kit (Ready to Use) (orb554687)

Description	Rat ATPase, Na ⁺ /K ⁺ Transporting Beta 4 Polypeptide (ATP1b4) Sandwich ELISA KitReady to Use: Pre-dilu...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9R193
Sensitivity	0.055ng/mL
Entrez	84396
Expiration Date	Please enquire.