

## Product Datasheet

### Rat bEP ELISA Kit (Ready to Use) (orb554696)

|                          |                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Rat Beta-Endorphin (bEP) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenie...                                                                        |
| <b>Reactivity</b>        | Rat                                                                                                                                                                            |
| <b>Range</b>             | 15.6-1000pg/mL                                                                                                                                                                 |
| <b>Form/Appearance</b>   | Ready to Use                                                                                                                                                                   |
| <b>Storage</b>           | All reagents should be stored at -20°C upon receipt                                                                                                                            |
| <b>Note</b>              | For research use only                                                                                                                                                          |
| <b>Application notes</b> | Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates or other biological fluids. |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernatants and other biological fluids.                                                   |
| <b>Sensitivity</b>       | 4.35pg/mL                                                                                                                                                                      |
| <b>Expiration Date</b>   | Please enquire.                                                                                                                                                                |