

Product Datasheet

Rat CASP8AP2 ELISA Kit (Ready to Use) (orb554754)

Description	Rat Caspase 8 Associated Protein 2 (CASP8AP2) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D4A7V6
Sensitivity	0.136ng/mL
Entrez	313128
Expiration Date	Please enquire.