

Product Datasheet

Rat CR2 ELISA Kit (Ready to Use) (orb554822)

Description	Rat Complement Receptor 2 (CR2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	A0A0G2K684
Sensitivity	0.57ng/mL
Entrez	289395
Expiration Date	Please enquire.