

Product Datasheet

Rat DIDO1 ELISA Kit (Ready to Use) (orb554859)

Description	Rat Death Inducer Obliterator 1 (DIDO1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	D3ZWL9
Sensitivity	0.059ng/mL
Entrez	362286
Expiration Date	Please enquire.