

Product Datasheet

Rat DEFb2 ELISA Kit (Ready to Use) (orb554862)

Description	Rat Defensin Beta 2 (DEFb2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conve...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Sensitivity	0.061ng/mL
Expiration Date	Please enquire.