

Product Datasheet

Rat Dopamine D2 receptor ELISA Kit (Ready to Use) (orb554879)

Description	Rat Dopamine Receptor D2 (Dopamine D2 receptor) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at - 20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P61169
Sensitivity	0.117ng/mL
Entrez	24318
Expiration Date	Please enquire.