

Product Datasheet

Rat DUOX2 ELISA Kit (Ready to Use) (orb554884)

Description	Rat Dual Oxidase 2 (DUOX2) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9ES45
Sensitivity	0.058ng/mL
Entrez	79107
Expiration Date	Please enquire.