

Product Datasheet

Rat EML2 ELISA Kit (Ready to Use) (orb554886)

Description	Rat Echinoderm Microtubule Associated Protein Like Protein 2 (EML2) Sandwich ELISA KitReady to Use: ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6P6T4
Sensitivity	0.126ng/mL
Entrez	192360
Expiration Date	Please enquire.