

Product Datasheet

Rat EG-VEGF ELISA Kit (Ready to Use) (orb554893)

Description	Rat Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) Sandwich ELISA KitReady to ...
Reactivity	Rat
Range	31.2-2000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q8R414
Sensitivity	13.6pg/mL
Entrez	192205
Expiration Date	Please enquire.