

Product Datasheet

Rat NNE ELISA Kit (Ready to Use) (orb554901)

Description	Rat Enolase, Non Neuronal (NNE) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for c...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04764
Sensitivity	0.114ng/mL
Entrez	24333
Expiration Date	Please enquire.