

Product Datasheet

Rat ERa ELISA Kit (Ready to Use) (orb554912)

Description	Rat Estrogen Receptor Alpha (ERa) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P06211
Sensitivity	0.117ng/mL
Entrez	24890
Expiration Date	Please enquire.