

Product Datasheet

Rat EAAT2 ELISA Kit (Ready to Use) (orb554915)

Description	Rat Excitatory Amino Acid Transporter 2 (EAAT2) Sandwich ELISA KitReady to Use: Pre-diluted detectio...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P31596
Sensitivity	0.058ng/mL
Entrez	29482
Expiration Date	Please enquire.