

Product Datasheet

Rat FGF12 ELISA Kit (Ready to Use) (orb554948)

Description	Rat Fibroblast Growth Factor 12 (FGF12) Sandwich ELISA KitReady to Use: Pre-diluted detection reagen...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P61150
Sensitivity	5.9pg/mL
Entrez	170630
Expiration Date	Please enquire.