

Product Datasheet

CCR2 antibody (orb55496)

Description	Rabbit polyclonal antibody against CCR2 conjugated to Biotin.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	IHC, WB
Immunogen	Synthetic peptide of Human CCR2
Target	CCR2
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-704 was sufficient for detection of CCR2 in 15 µg of rat brain cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	P41597
NCBI	NP_001116513.2
Entrez	729230
Dilution Range	WB (1:1000); IHC (1:50)
Expiration Date	12 months from date of receipt.