

Product Datasheet

Rat gp130 ELISA Kit (Ready to Use) (orb554987)

Description	Rat Glucoprotein 130 (gp130) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conv...
Reactivity	Rat
Range	3.12-200ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P40190
Sensitivity	1.21ng/mL
Entrez	25205
Expiration Date	Please enquire.