


Product Datasheet

Insulin Receptor Beta (phospho-Tyr1185) antibody (orb5550)

Description	Insulin Receptor Beta (phospho-Tyr1185) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor Beta around the phosphorylation site of Tyr1185 (DI(p-Y)ET)
Target	Phospho-Insulin Receptor Beta (Tyr1185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68/152 kDa
CAS Number	3
Uniprot ID	P06213
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...