

Product Datasheet

Rat GSTa3 ELISA Kit (Ready to Use) (orb555013)

Description	Rat Glutathione S Transferase Alpha 3 (GSTa3) Sandwich ELISA KitReady to Use: Pre-diluted detection ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P46418
Sensitivity	0.079ng/mL
Entrez	494500
Expiration Date	Please enquire.