

Product Datasheet

Rat GDF2 ELISA Kit (Ready to Use) (orb555036)

Description	Rat Growth Differentiation Factor 2 (GDF2) Sandwich ELISA KitReady to Use: Pre-diluted detection rea...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	M0RB87
Sensitivity	0.126ng/mL
Entrez	290921
Expiration Date	Please enquire.