

Product Datasheet

Rat HSPb6 ELISA Kit (Ready to Use) (orb555054)

Description	Rat Heat Shock Protein Beta 6 (HSPb6) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P97541
Sensitivity	0.068ng/mL
Entrez	192245
Expiration Date	Please enquire.