

Product Datasheet

Rat Hepcidin-25 ELISA Kit (Ready to Use) (orb555067)

Description	Rat Hepcidin (Hepcidin-25) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for conven...
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q99MH3
Sensitivity	0.227ng/mL
Entrez	84604
Expiration Date	Please enquire.