

Product Datasheet

Rat INHA ELISA Kit (Ready to Use) (orb555097)

Description	Rat Inhibin A (INHA) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient a...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P17490
Sensitivity	6.5pg/mL
Entrez	24504
Expiration Date	Please enquire.