

Product Datasheet

Rat IL18BP ELISA Kit (Ready to Use) (orb555141)

Description	Rat Interleukin 18 Binding Protein (IL18BP) Sandwich ELISA KitReady to Use: Pre-diluted detection re...
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F7F5P4
Sensitivity	0.215ng/mL
Entrez	84388
Expiration Date	Please enquire.