

Product Datasheet

Rat IL2Ra ELISA Kit (Ready to Use) (orb555144)

Description	Rat Interleukin 2 Receptor Alpha (IL2Ra) Sandwich ELISA KitReady to Use: Pre-diluted detection reage...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P26897
Sensitivity	0.051ng/mL
Entrez	25704
Expiration Date	Please enquire.