

Product Datasheet

Rat Leptin Receptor ELISA Kit (Ready to Use) (orb555193)

Description	Rat Leptin Receptor (Leptin Receptor) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents...
Reactivity	Rat
Range	1.56-100ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q62959
Sensitivity	0.66ng/mL
Entrez	24536
Expiration Date	Please enquire.