

Product Datasheet

KSR1 (phospho-Ser392) antibody (orb5553)

Description	KSR1 (phospho-Ser392) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human KSR1 around the phosphorylation site of Ser392 (RI(p-S)FL)
Target	Phospho-KSR1 (Ser392)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102 kDa
Uniprot ID	Q8IVT5
Expiration Date	12 months from date of receipt.