

Product Datasheet

GNDF antibody (orb55536)

Description	Rabbit polyclonal antibody against GDNF conjugated to PerCP.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PerCP
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide from the mid-protein of Human GDNF
Target	GDNF
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-710 was sufficient for detection of GDNF in 15 µg of human HeLa cell lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	24 kDa
Uniprot ID	P39905
NCBI	NP_000505.1
Entrez	2668
Dilution Range	WB (1:1000), ICC/IF (1:100); IHC (1:50)
Expiration Date	12 months from date of receipt.