

Product Datasheet

Rat PRF1 ELISA Kit (Ready to Use) (orb555368)

Description	Rat Perforin 1 (PRF1) Sandwich ELISA Kit Ready to Use: Pre-diluted detection reagents for convenient ...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stability serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q5FVS5
Sensitivity	0.18ng/mL
Entrez	50669
Expiration Date	Please enquire.