

Product Datasheet

Rat RAC1 ELISA Kit (Ready to Use) (orb555460)

Description	Rat Ras Related C3 Botulinum Toxin Substrate 1 (RAC1) Sandwich ELISA KitReady to Use: Pre-diluted de...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q6RUV5
Sensitivity	0.069ng/mL
Entrez	363875
Expiration Date	Please enquire.