

Product Datasheet

Rat ARHGEF7 ELISA Kit (Ready to Use) (orb555479)

Description	Rat Rho Guanine Nucleotide Exchange Factor 7 (ARHGEF7) Sandwich ELISA KitReady to Use: Pre-diluted d...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	055043
Sensitivity	0.064ng/mL
Entrez	114559
Expiration Date	Please enquire.