

Product Datasheet

Rat SOX2 ELISA Kit (Ready to Use) (orb555510)

Description	Rat Sex Determining Region Y Box Protein 2 (SOX2) Sandwich ELISA KitReady to Use: Pre-diluted detect...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates or other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	D4A543
Sensitivity	0.064ng/mL
Entrez	499593
Expiration Date	Please enquire.