

Product Datasheet

Rat STAT3 ELISA Kit (Ready to Use) (orb555515)

Description	Rat Signal Transducer And Activator Of Transcription 3 (STAT3) Sandwich ELISA KitReady to Use: Pre-d...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P52631
Sensitivity	0.054ng/mL
Entrez	25125
Expiration Date	Please enquire.