

Product Datasheet

Rat SPB ELISA Kit (Ready to Use) (orb555560)

Description	Rat Surfactant Associated Protein B (SPB) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates, cell lysates, lung lavage fluid or other biological fluids.
Sample Types	tissue homogenates, cell lysates, lung lavage fluid and other biological fluids.
Uniprot ID	P22355
Sensitivity	0.142ng/mL
Entrez	192155
Expiration Date	Please enquire.