

Product Datasheet

Rat TAUT ELISA Kit (Ready to Use) (orb555571)

Description	Rat Taurine Transporter (TAUT) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for co...
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilitytissue homogenates or other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P31643
Sensitivity	0.052ng/mL
Entrez	29464
Expiration Date	Please enquire.