

Product Datasheet

Rat TNFRSF5 ELISA Kit (Ready to Use) (orb555631)

Description	Rat Tumor Necrosis Factor Receptor Superfamily, Member 5 (TNFRSF5) Sandwich ELISA KitReady to Use: P...
Reactivity	Rat
Range	15.6-1000pg/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma or other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q4QQW2
Sensitivity	5.2pg/mL
Entrez	171369
Expiration Date	Please enquire.