

Product Datasheet

Sheep IGF1 ELISA Kit (Ready to Use) (orb555685)

Description	Sheep Insulin Like Growth Factor 1 (IGF1) Sandwich ELISA KitReady to Use: Pre-diluted detection reagent...
Reactivity	Sheep
Range	0.234-15ng/mL
Form/Appearance	Ready to Use
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Application notes	Pre-diluted detection reagents and enhanced stabilityserum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P10763
Sensitivity	0.115ng/mL
Entrez	443318
Expiration Date	Please enquire.