

Product Datasheet

Primate PD-L2 protein (orb555718)

Description	Recombinant primate Programmed cell death ligand 2/PD-L2 (C-Fc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Leu20-Pro219
Source	Human Cells
Uniprot ID	A4GW30
Expiration Date	6 months from date of receipt.