

Product Datasheet

Primate CD47 protein (orb555739)

Description	Recombinant primate CD47 molecule/CD47 (C-Fc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50 mM Tris, 100 mM Glycine, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50 mM Tris, 100 mM Glycine, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Gln19-Pro139
Source	Human Cells
Uniprot ID	F7A802
Expiration Date	6 months from date of receipt.