

Product Datasheet

Primate CRTAM protein (orb555742)

Description	Recombinant primate Cytotoxic and regulatory T-cell molecule/CRTAM (C-Fc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 50 mM Tris, 100 mM Glycine, pH 7.5.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50 mM Tris-HCl, 100 mM Glycine, pH 7.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Ser18-Gly287
Source	Human Cells
Uniprot ID	A0A2K5TKL4
Expiration Date	6 months from date of receipt.