

Product Datasheet

Human Reg3protein (orb555744)

Description	Recombinant human Regenerating Islet-Derived Protein 3-γ/Reg3γ/REG3G (C-Fc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of PBS,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS,pH7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-Fc
Note	For research use only
Protein Sequence	Glu27-Asp175
Source	Human Cells
Uniprot ID	Q6UW15
Expiration Date	6 months from date of receipt.