

Product Datasheet

Human B-CK protein (orb555753)

Description	Recombinant human Creatine kinase B chain/B-CK (N-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,10% Glycerol,pH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl,150mM NaCl,10% Glycerol,pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Protein Sequence	Met1-Lys381
Source	E.coli
Uniprot ID	P12277
Expiration Date	6 months from date of receipt.