

Product Datasheet

Human VISTA protein (orb555770)

Description	Recombinant human Platelet Receptor Gi24/VISTA/B7-H5/PD-1H (C-mFc)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-mFc
Note	For research use only
Protein Sequence	Phe33-Ala194
Source	Human Cells
Uniprot ID	AAH20568.1
Expiration Date	6 months from date of receipt.