

Product Datasheet

Human TNFRSF10B protein (orb555801)

Description	Recombinant human TRAIL R2/TNFRSF10B/DR5/CD262 (C-Fc-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at - 20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc&His
Note	For research use only
Protein Sequence	Ile56-Glu182
Source	Human Cells
Uniprot ID	O14763
Expiration Date	6 months from date of receipt.